

STERIS: Certificate Of Processing

Prepared    STIRILAB  
For  
Ethylene Oxide Run ID 2002-4342

| Product Code       | Lot Number | Quantity (PL) |
|--------------------|------------|---------------|
| STIRILAB 03 EQ-15% | 202405075  | 18            |

Product meets Customer specifications; zero nonconformities occurred during this run.

Other Information

Sterilization Batch Number 22241102

Customer Specification ID: 37

**Signature Manifest**

Reviewed and E-Signed By: **Sarka Mejzlikova (Quality Analyst)**

DateTime Signed 04-Nov-2024 9:43 AM

Document Content Revision 1

Operating facilities are in compliance with applicable state and federal regulations (FDA, NRC, EPA, and OSHA) and provide services under a quality system which meets the requirements of FDA QSR, EN/ISO 13485, and in alignment with EN ANSI/AAMI/ISO 11135.

This Certificate becomes valid only with the Process Run Record Summary and Biological Indicators Incubation Certificate referring to applicable sterilization Batch Number.

Processing Location

STERIS AST CZ s.r.o.  
Průmyslová zóna Košíkov 80  
Velká Bíteš  
595 01

\*\*\* End of Ethylene Oxide Certificate of Processing Document \*\*\*

# Certificate of Analysis



Customer: Stiriab s.r.o.  
Address: Textilná 23  
Ružomberok  
SVK  
03401

Job No: BI-CZE-VEL-202402749  
Report Number: 130497  
Date Issued: 11-Nov-2024  
Page 1 of 1

Purchase Order No: 202405075  
Product: STIRILAB 03 EQ-15%, Biological Indicator  
Description: SGM Spore Strip – Bacillus atrophaeus ATCC #9372  
Standard Method: ISO 11138-1:2017  
Internal Work Instruction: CZE VEL : LWM-LAB002  
Customer Protocol: STIRILAB 03 EQ-15%

Date Sample(s) Received: 04-Nov-2024  
Date On Test: 04-Nov-2024

## Analytical Results

Label ID / Designation: 202405075

Lot / Sterilization Batch: 22241102

Microbiology

Result

Biological Indicators

|                          |                    |
|--------------------------|--------------------|
| BI Lot Number            | 488                |
| Number of samples        | 28                 |
| Positive Control Check 7 | Positive           |
| Number of positives      | 0                  |
| Negative Controls        | Negative           |
| Time on test             | 04/NOV/24 10:15:05 |
| Time off Test            | 11/NOV/24 10:23:57 |

The results presented relate only to the items tested above as supplied by the Customer and received by the laboratory.

This test report should not be reproduced except in full, without prior written permission.

Released by:

Eva Tomkova

Laboratory Operation Manager

Authorisation Date: 11-November-2024

Testing Facility  
CZE Velka Bites  
STERIS AST CZ s.r.o.  
Prumyslová zóna Košíkov 80

Velká Bíteš

595 01

# Process Run Record

## Summary

**STERIS AST CZ s.r.o.**  
Průmyslová zóna Košíkov 80  
595 01 Velká Bíteš, Czech Republic  
Tel: +420566562145  
Fax:

|                |                            |               |        |            |                    |
|----------------|----------------------------|---------------|--------|------------|--------------------|
| Run ID:        | 317157                     | Pallet Count: | 18     | Batch No.: | 22241102           |
| Customer:      | routine_STIRILAB 03 EQ-15% | Checksum:     | -14460 | Method:    | STIRILAB 03 EQ-15% |
| Specification: | 1.4                        |               |        |            |                    |

## 1.Preconditioning

|               |                    |                 |                     |
|---------------|--------------------|-----------------|---------------------|
| Equipment:    | VB Precond G       | Config. Rev.:   | Pending             |
| Date Started: | 01.11.2024 7:47:33 | Date Completed: | 01.11.2024 19:07:06 |
|               |                    | Duration:       | 11:19:33            |

### 1.1.Preconditioning

Phase No. 1

| Device           | Summary      | Summary Value | Specified Range        | ✓/✗ |
|------------------|--------------|---------------|------------------------|-----|
| Total Phase Time | Final - P104 | 11:19:33      | 11:00:00 - 1d 00:00:00 | ✓   |
| Time-InSpec      | Final - P104 | 11:05:00      | 11:00:00 - 1d 00:00:00 | ✓   |
| TE-Min-In-Spec   | Min - P101   | 47,4 °C       | 45,0 - 55,0 °C         | ✓   |
| TempAve          | Max - P101   | 51,3 °C       | 45,0 - 55,0 °C         | ✓   |
| RH-Min-In-Spec   | Min - P22    | 55,1 %        | 45,0 - 95,0 %          | ✓   |
| AT-RH-Avg        | Max - P22    | 64,0 %        | 45,0 - 95,0 %          | ✓   |

### 1.2.Transfer to Vessel

Phase No. 2

| Device            | Summary      | Summary Value | Specified Range | ✓/✗ |
|-------------------|--------------|---------------|-----------------|-----|
| Transfer Duration | Final - P104 | 00:17:07      | n/a - 00:45:00  | ✓   |

### Stage Summary

| Steps | Deviations | Date Started       | Date Ended          | Total Duration: |
|-------|------------|--------------------|---------------------|-----------------|
| 2     | 0          | 01.11.2024 7:47:33 | 01.11.2024 19:24:13 | 11:19:33        |

## Process Run Record Summary

Run ID 317157

## 2.Vessel Processing

Equipment: VB Chamber 22

Config. Rev.: 4

Date Started: 01.11.2024 19:24:13

Date Completed: 02.11.2024 6:49:58

Duration: 11:25:45

### 2.1.Start-up Phase

Phase No. 1

| Device   | Summary    | Summary Value | Specified Range | ✓/✗ |
|----------|------------|---------------|-----------------|-----|
| TE-CH-Hi | Final - P3 | 49,9 °C       | 48,0 - 54,0 °C  | ✓   |
| TE-CH-Lo | Final - P3 | 49,9 °C       | 48,0 - 54,0 °C  | ✓   |

### 2.2.Vacuum

Phase No. 2

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 00:13:24      | 00:10:00 - 00:45:00 | ✓   |
| Total Phase Time   | Final - P11 | 00:13:24      | 00:10:00 - 00:45:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 49,5 °C       | 40,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 49,9 °C       | 40,0 - 55,0 °C      | ✓   |
| PCurr              | Final - P14 | 22,0 kPa      | 12,0 - 32,0 kPa     | ✓   |

### 2.3.Leak Check

Phase No. 3

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 00:05:01      | 00:01:00 - 00:20:00 | ✓   |
| Total Phase Time   | Final - P11 | 00:05:01      | 00:01:00 - 00:20:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 49,8 °C       | 40,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 50,2 °C       | 40,0 - 55,0 °C      | ✓   |
| PCurr              | Min - P14   | 22,0 kPa      | 12,0 - 33,0 kPa     | ✓   |
| PCurr              | Max - P14   | 22,6 kPa      | 12,0 - 33,0 kPa     | ✓   |

### 2.4.Vacuum

Phase No. 4

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 00:00:21      | 00:00:01 - 00:05:00 | ✓   |
| Total Phase Time   | Final - P11 | 00:00:21      | 00:00:01 - 00:05:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 50,2 °C       | 40,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 50,2 °C       | 40,0 - 55,0 °C      | ✓   |
| PCurr              | Final - P14 | 22,0 kPa      | 21,0 - 22,2 kPa     | ✓   |

### 2.5.Steam Inject by Pressure

Phase No. 5

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 00:03:52      | 00:01:00 - 00:12:00 | ✓   |
| Total Phase Time   | Final - P11 | 00:03:52      | 00:01:00 - 00:12:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 50,2 °C       | 45,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 50,4 °C       | 45,0 - 55,0 °C      | ✓   |
| PCurr              | Final - P14 | 23,1 kPa      | 23,0 - 26,6 kPa     | ✓   |



## Process Run Record Summary

Run ID 317157

**Stage »** Vessel Processing

### 2.6. Conditioning Dwell by Pressure

Phase No. 6

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 01:25:01      | 01:10:00 - 01:40:00 | ✓   |
| Total Phase Time   | Final - P11 | 01:25:01      | 01:10:00 - 01:40:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 50,0 °C       | 45,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 50,7 °C       | 45,0 - 55,0 °C      | ✓   |
| PCurr              | Min - P14   | 23,0 kPa      | 22,0 - 26,6 kPa     | ✓   |
| PCurr              | Max - P14   | 23,5 kPa      | 22,0 - 26,6 kPa     | ✓   |
| RH Reported        | Min - P22   | 61,1 %        | N/A - N/A %         |     |
| RH Reported        | Max - P22   | 82,8 %        | N/A - N/A %         |     |

### 2.7. Vacuum

Phase No. 7

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 00:00:38      | 00:00:01 - 00:05:00 | ✓   |
| Total Phase Time   | Final - P11 | 00:00:38      | 00:00:01 - 00:05:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 49,9 °C       | 40,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 50,0 °C       | 40,0 - 55,0 °C      | ✓   |
| PCurr              | Final - P14 | 22,0 kPa      | 21,0 - 22,2 kPa     | ✓   |

### 2.8. EO Inject by Pressure

Phase No. 8

| Device              | Summary       | Summary Value | Specified Range     | ✓/✗ |
|---------------------|---------------|---------------|---------------------|-----|
| In-Spec Phase Time  | Final - P11   | 01:04:17      | 00:45:00 - 01:30:00 | ✓   |
| Total Phase Time    | Final - P11   | 01:04:17      | 00:45:00 - 01:30:00 | ✓   |
| TE-CH-Lo            | Min - P3      | 49,9 °C       | 46,0 - 54,0 °C      | ✓   |
| TE-CH-Hi            | Max - P3      | 50,8 °C       | 46,0 - 54,0 °C      | ✓   |
| PCurr               | Final - P1028 | 215,0 kPa     | n/a - N/A kPa       |     |
| PCurr               | Final - P14   | 215,0 kPa     | 200,0 - 225,0 kPa   | ✓   |
| EO Conc. Reported   | Final - P33   | 429,3 mg/l    | N/A - N/A mg/l      |     |
| EO Conc. Calculated | Final - P33   | 475,2 mg/l    | N/A - N/A mg/l      |     |

**Process Run Record Summary**

Run ID 317157

**Stage » Vessel Processing****2.9.Super-Atmospheric EO Dwell**

Phase No. 9

| Device              | Summary     | Summary Value | Specified Range     | ✓/✗ |
|---------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time  | Final - P11 | 06:00:01      | 05:55:00 - 07:00:00 | ✓   |
| Total Phase Time    | Final - P11 | 06:00:01      | 05:55:00 - 07:00:00 | ✓   |
| TE-CH-Lo            | Min - P3    | 49,4 °C       | 46,0 - 54,0 °C      | ✓   |
| TE-CH-Hi            | Max - P3    | 50,3 °C       | 46,0 - 54,0 °C      | ✓   |
| PCurr               | Min - P14   | 210,3 kPa     | 200,0 - 225,0 kPa   | ✓   |
| PCurr               | Min - P1028 | 210,3 kPa     | N/A - N/A kPa       |     |
| PCurr               | Max - P1028 | 215,0 kPa     | N/A - N/A kPa       |     |
| PCurr               | Max - P14   | 215,0 kPa     | 200,0 - 225,0 kPa   | ✓   |
| RH Reported         | Min - P22   | 0,0 %         | N/A - N/A %         |     |
| RH Reported         | Max - P22   | 0,0 %         | N/A - N/A %         |     |
| EO Conc. Reported   | Min - P33   | 304,4 mg/l    | N/A - N/A mg/l      |     |
| EO Conc. Reported   | Max - P33   | 429,3 mg/l    | N/A - N/A mg/l      |     |
| EO Conc. Calculated | Min - P33   | 475,0 mg/l    | N/A - N/A mg/l      |     |
| EO Conc. Calculated | Max - P33   | 476,3 mg/l    | N/A - N/A mg/l      |     |

**2.10.Vacuum**

Phase No. 10

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 00:22:55      | 00:15:00 - 00:45:00 | ✓   |
| Total Phase Time   | Final - P11 | 00:22:55      | 00:15:00 - 00:45:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 48,1 °C       | 40,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 50,1 °C       | 40,0 - 55,0 °C      | ✓   |
| PCurr              | Final - P14 | 26,0 kPa      | 16,0 - 36,0 kPa     | ✓   |

**2.11.Air Inbleed**

Phase No. 11

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 00:07:21      | 00:03:00 - 00:17:00 | ✓   |
| Total Phase Time   | Final - P11 | 00:07:21      | 00:03:00 - 00:17:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 48,2 °C       | 40,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 48,6 °C       | 40,0 - 55,0 °C      | ✓   |
| PCurr              | Final - P14 | 93,0 kPa      | 90,0 - 103,0 kPa    | ✓   |

**2.12.Vacuum**

Phase No. 12

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 00:12:08      | 00:12:00 - 00:28:00 | ✓   |
| Total Phase Time   | Final - P11 | 00:12:08      | 00:12:00 - 00:28:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 47,8 °C       | 40,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 48,6 °C       | 40,0 - 55,0 °C      | ✓   |
| PCurr              | Final - P14 | 27,0 kPa      | 17,0 - 37,0 kPa     | ✓   |

**Process Run Record Summary**

Run ID 317157

**Stage » Vessel Processing****2.13. Air Inbleed**

Phase No. 13

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 00:07:13      | 00:03:00 - 00:17:00 | ✓   |
| Total Phase Time   | Final - P11 | 00:07:13      | 00:03:00 - 00:17:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 48,1 °C       | 40,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 49,6 °C       | 40,0 - 55,0 °C      | ✓   |
| PCurr              | Final - P14 | 93,0 kPa      | 90,0 - 103,0 kPa    | ✓   |

**2.14. Vacuum**

Phase No. 14

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 00:12:08      | 00:12:00 - 00:28:00 | ✓   |
| Total Phase Time   | Final - P11 | 00:12:08      | 00:12:00 - 00:28:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 48,4 °C       | 40,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 49,6 °C       | 40,0 - 55,0 °C      | ✓   |
| PCurr              | Final - P14 | 27,0 kPa      | 17,0 - 37,0 kPa     | ✓   |

**2.15. Air Inbleed**

Phase No. 15

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 00:07:14      | 00:03:00 - 00:17:00 | ✓   |
| Total Phase Time   | Final - P11 | 00:07:14      | 00:03:00 - 00:17:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 48,6 °C       | 40,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 50,6 °C       | 40,0 - 55,0 °C      | ✓   |
| PCurr              | Final - P14 | 93,0 kPa      | 90,0 - 103,0 kPa    | ✓   |

**2.16. Vacuum**

Phase No. 16

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 00:12:08      | 00:12:00 - 00:28:00 | ✓   |
| Total Phase Time   | Final - P11 | 00:12:08      | 00:12:00 - 00:28:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 49,0 °C       | 40,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 50,6 °C       | 40,0 - 55,0 °C      | ✓   |
| PCurr              | Final - P14 | 27,0 kPa      | 17,0 - 37,0 kPa     | ✓   |

**2.17. Air Inbleed**

Phase No. 17

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 00:07:15      | 00:03:00 - 00:17:00 | ✓   |
| Total Phase Time   | Final - P11 | 00:07:15      | 00:03:00 - 00:17:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 49,1 °C       | 40,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 51,0 °C       | 40,0 - 55,0 °C      | ✓   |
| PCurr              | Final - P14 | 93,1 kPa      | 90,0 - 103,0 kPa    | ✓   |

**Process Run Record      Summary**

Run ID      317157

**Stage »    Vessel Processing****2.18.End Phase**

Phase No. 18

| Device             | Summary     | Summary Value | Specified Range     | ✓/✗ |
|--------------------|-------------|---------------|---------------------|-----|
| In-Spec Phase Time | Final - P11 | 01:04:43      | 00:00:01 - 05:00:00 | ✓   |
| Total Phase Time   | Final - P11 | 01:04:43      | 00:00:01 - 05:00:00 | ✓   |
| TE-CH-Lo           | Min - P3    | 50,0 °C       | 40,0 - 55,0 °C      | ✓   |
| TE-CH-Hi           | Max - P3    | 51,0 °C       | 40,0 - 55,0 °C      | ✓   |
| PCurr              | Min - P14   | 93,0 kPa      | 90,0 - 103,0 kPa    | ✓   |
| PCurr              | Max - P14   | 94,4 kPa      | 90,0 - 103,0 kPa    | ✓   |

**Stage Summary**

| Steps | Deviations | Date Started        | Date Ended         | Total Duration: |
|-------|------------|---------------------|--------------------|-----------------|
| 18    | 0          | 01.11.2024 19:24:13 | 02.11.2024 6:49:58 | 11:25:45        |



3.Aeration

|               |                    |                 |                     |
|---------------|--------------------|-----------------|---------------------|
| Equipment:    | VB Aeration 04     | Config. Rev.:   | Pending             |
| Date Started: | 02.11.2024 7:18:59 | Date Completed: | 02.11.2024 16:53:33 |
|               |                    | Duration:       | 09:34:34            |

3.1.Aeration

Phase No. 1

| Device           | Summary      | Summary Value | Specified Range        | ✓/✗ |
|------------------|--------------|---------------|------------------------|-----|
| Total Phase Time | Final - P104 | 09:34:34      | 00:00:00 - 7d 00:00:00 | ✓   |
| Time-InSpec      | Final - P104 | 09:30:00      | 00:00:00 - 7d 00:00:00 | ✓   |
| TE-Min-In-Spec   | Min - P101   | 45,7 °C       | 41,0 - 51,0 °C         | ✓   |
| TempAve          | Max - P101   | 46,5 °C       | 41,0 - 51,0 °C         | ✓   |

Stage Summary

| Steps | Deviations | Date Started       | Date Ended          | Total Duration: |
|-------|------------|--------------------|---------------------|-----------------|
| 1     | 0          | 02.11.2024 7:18:59 | 02.11.2024 16:53:33 | 09:34:34        |

Minimum Load Temperature before PRECOND: ≥ 10,8 ° C value: 21,2 ° C  
SIADTOX 15% - specification: N/A value: 119,9 kg  
CYLINDER NR.  
1. 8422  
2. 161  
3. 3213  
4. 2739

CONFORMING

Vrbka,  
Vojtech

Digitálně podepsal  
Vrbka, Vojtech  
Datum: 2024.11.04  
08:07:19 +01'00'



## Transport Document



202405075

STIRILAB s.r.o. a dominique Dutscher company

Textilná 23

034 01

Ružomberok

VAT No. SK2021729545

Tel +421 44 432 27 20 Fax +421 44 432 27 17

www.stirilab.com, email ekonom@stirilab.com

IBAN : SK 311111 000000 102498 0028, SWIFT : UNCRSKBX

## CUSTOMER

STERIS AST CZ s.r.o.

Průmyslová zóna Košíkov 80

595 01 Velká Bíteš, ČR

## SHIP TO

STERIS AST CZ s.r.o.

Průmyslová zóna Košíkov 80

59501 Velká Bíteš Česká republika

| CUST. NO. | VAT CODE   | DATE       | TD NUMBER | ISSUED BY | PAGE   |
|-----------|------------|------------|-----------|-----------|--------|
|           | CZ27159621 | 31.10.2024 | 202405075 |           | 1 to 2 |

| CODE    | DESCRIPTION                              | QTY per BOX    | MU        | QUANTITY |
|---------|------------------------------------------|----------------|-----------|----------|
| BSC740  | Fertility collector                      | 40             | 40,000    | 41,00    |
| order   | lotto 24P2004309                         | exp 31.10.2028 |           |          |
| BSC740  | Fertility collector                      | 40             | 40,000    | 23,00    |
| order   | lotto 24P2004310                         | exp 31.10.2028 |           |          |
| BSV120  | ANSA 10 uL, PP, peel pack 20 pcs,sterile | 10000          | 1 000,000 | 70,00    |
| order   | lotto 24P2003407                         | exp 31.08.2028 |           |          |
| BSV120  | ANSA 10 uL, PP, peel pack 20 pcs,sterile | 10000          | 1 000,000 | 740,00   |
| order   | lotto 24P2003409                         | exp 31.10.2028 |           |          |
| BSV1200 | ANSA 10 uL, PP, peel pack 10 pcs,sterile | 12000          | 1 000,000 | 72,00    |
| order   | lotto 24P2003304                         | exp 31.10.2028 |           |          |
| BSV1201 | ANSA 10 uL, sterile                      | 2100           | 2 100,000 | 32,00    |
| order   | lotto 24P2003707                         | exp 31.10.2028 |           |          |
| BSV1205 | Anse jednorázové, 10ul                   | 4000           | 1 000,000 | 16,00    |
| order   | lotto 24P2003205                         | exp 31.10.2028 |           |          |
| BSV121  | ANSA 1 uL, PP, peel pack 20 pcs,sterile  | 10000          | 1 000,000 | 800,00   |
| order   | lotto 24P2003010                         | exp 31.10.2028 |           |          |
| BSV121  | ANSA 1 uL, PP, peel pack 20 pcs,sterile  | 10000          | 1 000,000 | 680,00   |
| order   | lotto 24P2003010                         | exp 31.10.2028 |           |          |
| BSV1210 | Anse single-use, volume 1 ul             | 12000          | 1 000,000 | 24,00    |
| order   | lotto 24P2003502                         | exp 31.10.2028 |           |          |
| BSV1215 | ANSA 1 uL, PP, peel pack 5 pcs, sterile  | 0              | 1 000,000 | 28,00    |
| order   | lotto 24P2003104                         | exp 31.10.2028 |           |          |



31.10.2024

19:00



**STIRILAB s.r.o.** *a dominique Dutscher company*

**Textilná 23**

**034 01**

**Ružomberok**

**VAT No. SK2021729545**

Tel +421 44 432 27 20 Fax +421 44 432 27 17

www.stirilab.com, email ekonom@stirilab.com

IBAN : SK 311111 000000 102498 0028, SWIFT : UNCRSKBX

## Transport Document

CUSTOMER

**STERIS AST CZ s.r.o.**

Průmyslová zóna Košíkov 80

595 01 Velká Bíteš, ČR

SHIP TO

**STERIS AST CZ s.r.o.**

Průmyslová zóna Košíkov 80

59501 Velká Bíteš Česká republika

| CUST. NO. | VAT CODE   | DATE       | TD NUMBER | ISSUED BY | PAGE   |
|-----------|------------|------------|-----------|-----------|--------|
|           | CZ27159621 | 31.10.2024 | 202405075 |           | 2 to 2 |

| CODE | DESCRIPTION | QTY per BOX | MU | QUANTITY |
|------|-------------|-------------|----|----------|
|------|-------------|-------------|----|----------|

SHIPPING BY

EXW

ISSUED BY:

Božeková

TRANSPORT

SALE

TRANSPORT START DATE

31.10.2024

DELIVERY TERM

EXW

BOXES

344

PALLETS

18

GROSS WEIGHT KG

2 362,10

NETTO WEIGHT

2 140,99

VOLUME m3

18,14

DRIVER'S SIGNATURE

ANNOTATIONS

RECIPIENT'S SIGNATURE

*Any complaint must be made within 8 days after goods receipt. Goods are considered of customer property only after full payment receipt. Any disputes between the are subject to the exclusive jurisdiction of the courts of Slovak republic*

**ACCEPTANCE  
VERIFICATION**